

Mouse SLA/LP protein, His tag

Cat. No. GTX57386-pro

Species Mouse**Package**

20 µg

Applications

Observed MW (kDa) 55 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 200mM NaCl, 50% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMPESFA AGERRVSPAY VRQGCEARRA
HEHLIRLLLE QGKCPEDGWD ESTLEFLHE LAVMDSNNFL GNCVGGEREG RVASALVARR HYRFIHGIGR SGDISAVQPK
AAGSSLLNKI TNSLVNLVIK LAGVHSVASC FVPMATGMS LTLCFLTLRH KRPKAKYIIW PRIDQKSCFK SMVTAGFEPV
VIENVLEGDE LRTDLKAVEA KIQELGPEHI LCLHSTTACF APRVPDRLEE LAVICANYDI PHVNNAYGL QSSKCMHLIQ
QGARVGRIDA FVQSLDKNFM VPVGGAIAG FNEPFIQDIS KMYPGRASAS PSLDVLITLL SLGCSGYRKL LKERKEMFVY
LSTQLKKLAE AHNERLLQTP HNPISLAMTL KTIDGHHDKA VTQLGSMLFT RQVSGARAVP LGNVQTVSGH TFRGFMSHAD
NYPCAYLNAA AAIGMKMQDV DLFIKRLDKC LNIVRKEQTR ASVWSGADR N KAEDADIEEM ALKLDDVLGD VGQGPAL

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57386-pro Image

3μg Mouse SLA/LP protein (GTX57386-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).