

Mouse CNDP1 protein, His tag

Cat. No. GTX57411-pro

Species Mouse

Package
50 µg

Applications

Observed MW (kDa) 56 kDa.

Properties

Form Liquid

Buffer PBS, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, C-terminal His-Tag: MFSSAHSGLL EKLPHYDLH QDEFVQTLKE WVAIESDSVQ PVPRLRQKLF QMMALAADKL RNLGAGVESI DLGSQQMPDG QSLPIPPILL AELGSDPEKP TVCFYGHLDV QPAQKDDGWL TDPYTLTEVD GKLYGRGATD NKGVLAWIN AVSTFRALQQ DLPVNIKLIL EGMEEAGSIA LEELVMREKD HFFSSVDYIV ISDNLWLSQR KPALTYGTRG NCYFTVEVKC RDQDFHSGTF GGILNEPMAD LVALLGSLVD SSGHILIPGI YDQMAPITEG EKTMYKNIDM DLEEQNINQ VEKFLFDKE ELLMHLWRYP SLSIHGIEGA FDEPGTKTVI PGRVLGKFSI RLVPTMSPSV VEKQVTQHLE AVFSKRNSFN KMAVSMVLGL HPWTANVNDT QYLAAQRTIK TVFGVNPDMI RDGSTIPIAK IFQAITQKSV MMLPLGAVDD GEHSQNEKIN RWNVIQGSKL FAAFFLELSK QHSGHQMPSS VYLEHHHHHHH

Expression System Baculovirus

Purity > 95% by SDS-PAGE.

Endotoxin < 1.0 EU/µg (determined by LAL method)

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES



GTX57411-pro Image

3μg Mouse CNDP1 protein (GTX57411-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).