

Human CUEDC1 protein, His tag

Cat. No. GTX57433-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 42 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMTSLFRR SSSGSGGGGT AGARGGGGGT
AAPQELNNSR PARQVRRLEF NQAMDDFKTM FPNMDYDIE CVLRANS GAV DATIDQLQM NLEGGGSSGG VYEDSSDSED
SIPPEILERT LEPDSSDEEP PPVYSPPAYH MHVFDRPYPL APPTPPPRID ALGSGAPTSQ RRYRNWNPPPL LGNLPDDFLR ILPQQQLDSIQ
GNAGGPKPGS GEGCPPAMAG PGPGDQESRW KQYLEDERIA LFLQNEEFMK ELQRNRDFLL ALERDRKYE SQKSKSSSVA
VGNDFGFSSP VPGTGDANPA VSEDALFRDK LKHMKGSTRR KLFELARAFS ETKMRKSKR KHLKHKQSLG AAASTANLLD
DVEGHACDED FRGRRQEAPK VEEGLREGQ

Expression System E. coli**Purity** > 85% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57433-pro Image

3μg Human CUEDC1 protein (GTX57433-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).