

Human PEF1 protein, His tag

Cat. No. GTX57443-pro

Species Human

Package

50 µg

Applications

Observed MW (kDa) 30 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 200mM NaCl, 50% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.25 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMASY PYRQGCPGAA
GQAPGAPPGS YYPGPPNSGG QYGSGLPPGG GYGGPAPGGP YGPPAGGGPY GHPNPGMFPS GTPGGPYGGA APGGPYGQPP
PSSYGAQQPG LYQGQGAPPN VDPEAYSWFQ SVDSDHSGYI SMKELKQALV NCNWSSFNDE TCLMMINMFD KTKSGRIDVY
GFSALWKFIQ QWKNLFQQYD RDRSGSISYT ELQQALSQMG YNLSPQFTQL LVSRYCPRSA NPAMQLDRFI QVCTQLQVLT
EAFREKDTAV QGNIRLSFED FVTMTASRML

Expression System E. coli

Purity > 85% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57443-pro Image

3μg Human PEF1 protein (GTX57443-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).