

Human NCF1 protein, His tag

Cat. No. GTX57451-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 45 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 100mM NaCl, 20% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, C-terminal His-Tag; MGDTFIRHIA LLGF EKRFVP SQHYVYMFLV KWQDLSEKW YRRFTEIYEF HKTLKEMFPI
EAGAINPENR IIPHLPAKWW FDGQRAAENR QGTLTEYCST LMSLPTKISR CPHLLDFFKV RPDDLKLPD NQTKKPETYL
MPKDGKSTAT DITGPILQT YRAIANYEKT SGSEMA LSTG DVVEVVEKSE SGWWFCQMKA KRGWIPASFL EPLDSPDETE
DPEPNYAGEP YVAIKAYTAV EGDEVSLLEG EAVEVIHKLL DGWWVIRKDD VTGYFPSMYL QKSGQDVSQA QRQIKRGAPP
RRSSIRNAHS IHQSRKRRLS QDAYRRNSVR FLQQRRRQAR PGPQSPGSPL EERQTQRSK PQPAVPPRPS ADLILNRCSE
STKRKLASAV VEH HHHHHH

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57451-pro Image

3μg Human NCF1 protein (GTX57451-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).