

## E. coli SecB protein

Cat. No. GTX57491-pro

**Species** E. coli**Package**  
100 µg

## Applications

**Observed MW (kDa)** 17 kDa.

## Properties


**Form** Liquid**Buffer** 20mM Tris-HCl, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full length protein, MSEQNNTMT FQIQRIYTKD ISFEAPNAPH VFQKDWQPEV KLDLDTASSQ LADDVYEWL RVTVTASLGE  
ETAFLCEVQQ GGIFSIAIE GTQMAHCLGA YCPNLFPPA RECITSMVSR GTFPQLNLAP VNFDALFMNY LQQQAGEGTE EHQDA**Expression System** E. coli**Purity** > 95% by SDS-PAGE.**Conjugation** Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



## GTX57491-pro Image

3µg SecB protein (GTX57491-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

For full product information, images and publications, please visit our [website](https://www.genetex.com).