

Human NANOGP8 protein, His tag

Cat. No. GTX57518-pro

Species Human

Package

20 µg

Applications

Observed MW (kDa) 35 kDa.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 0.4M Urea, 5% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, C-terminal His-Tag: MSVDPACQQS LPCFEASDCK ESSPMPVICG PEENYPQLQM SSAEMPHDET
VSPLPSSMDL LIQDSPDSST SPKGGKQPTSA ENSVAKKEDK VPKKKQKTRT VFSSTQLCVL NDRFQRQKYL SLQQMQELSN
ILNLSYKQVK TWFNQRMKS KRWQKNNWPK NSNGVTQKAS APTYPSLYSS YHGGCLVNPT GNLPMTWSNQT WNNSTWSNQT
QNIQSWSNHS WNTQTWCTQS WNNQAWNSPF YNCGEESLQS CMHFQPNSPA SDLEAALEAA GEGLNVIQQT TRYFSTPQTM
DLFLNYSMMN QPEDV

Expression System E. coli

Purity > 85% by SDS-PAGE.

Conjugation Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX57518-pro Image

3μg Human NANOGP8 protein (GTX57518-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).