

Pyruvate Kinase (liver/RBC) antibody [AT1E3]

Cat. No. GTX57581

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 62 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT1E3 is derived from hybridization of mouse F2 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PKLR protein.
Purification	Protein G Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57581 WB Image

WB analysis of various samples using GTX57581 Pyruvate Kinase (liver/RBC) antibody.

Lane 1 : Balb/3T3 whole cell lysate

Lane 2 : HeLa whole cell lysate


Lane 3 : Raji whole cell lysate

Lane 4 : K562 whole cell lysate

Lane 5 : A549 whole cell lysate

Loading : 40 µg

Dilution : 1:500




GTX57581 ICC/IF Image

ICC/IF analysis of Hep3B cells using GTX57581 Pyruvate Kinase (liver/RBC) antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100



GTX57581 WB Image

WB analysis of 293T whole cell lysate using GTX57581 Pyruvate Kinase (liver/RBC) antibody.

Loading : 40 µg

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).