

## PDCL3 antibody [AT8F9]

Cat. No. GTX57608

|              |            |
|--------------|------------|
| Host         | Mouse      |
| Clonality    | Monoclonal |
| Isotype      | IgG1       |
| Applications | WB, ICC/IF |
| Reactivity   | Human      |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution                     |
|--------------------|------------------------------------------|
| WB                 | Recommended starting dilution is 1:1000. |
| ICC/IF             | Assay dependent                          |

Not tested in other applications.

Calculated MW 28 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 10% Glycerol                                                                                                                                                                                                          |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | The clone AT8F9 is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PDCL5 protein.                                                               |
| Purification  | Protein A Purified                                                                                                                                                                                                         |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTx57608 WB Image

WB analysis of various samples using GTx57608 PDCL3 antibody.

Samples : Jurkat and A549 whole cell lysate

Loading : 40 µg

Dilution : 1:1000




### GTx57608 WB Image

WB analysis of Recombinant Human PDCL3 protein using GTx57608 PDCL3 antibody.

Loading : 100 ng

Dilution : 1:1000



### GTx57608 ICC/IF Image

ICC/IF analysis of HeLa cells using GTx57608 PDCL3 antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).