

ACAT2 antibody [AT55D5]

Cat. No. GTX57610

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Recommended starting dilution is 1:1000.
ICC/IF	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Calculated MW 41 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT55D5 is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ACAT4 protein.
Purification	Protein A Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX57610 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX57610 ACAT2 antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100


GTX57610 WB Image


WB analysis of various samples using GTX57610 ACAT2 antibody.

Lane 1 : 293T cell lysate

Lane 2 : ACAT2 transfected 293T cell lysate

Loading : 5 µg

Dilution : 1:1000


GTX57610 WB Image

WB analysis of various samples using GTX57610 ACAT2 antibody.

Lane 1 : HepG2 whole cell lysate

Lane 2 : Ramos whole cell lysate

Lane 3 : MCF-7 whole cell lysate


Lane 4 : 293T whole cell lysate

Lane 5 : HL-60 whole cell lysate

Lane 6 : SW480 whole cell lysate

Loading : 40 µg

Dilution : 1:1000


GTX57610 FCM Image

FACS analysis of PC-3 cells using GTX57610 ACAT2 antibody.

Cell Number: 1×10^6 cells

Upper panel : negative control, Lower panel: Primary antibody

Antibody amount: 2-5 µg



For full product information, images and publications, please visit our [website](http://www.genetex.com).