

Fumarate hydratase antibody [AT3F6]

Cat. No. GTX57625

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Calculated MW 55 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT3F6 is derived from hybridization of mouse F2 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human Fumarase protein.
Purification	Protein A Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX57625 ICC/IF Image

ICC/IF analysis of A549 cells using GTX57625 Fumarate hydratase antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100




GTX57625 FCM Image

FACS analysis of HeLa cells using GTX57625 Fumarate hydratase antibody.

Cell Number: 1×10^6 cells

Upper panel : negative control, Lower panel: Primary antibody

Antibody amount: 2-5 μ g



GTX57625 WB Image

WB analysis of various samples using GTX57625 Fumarate hydratase antibody.

Lane 1 : HeLa whole cell lysate

Lane 2 : LNCap whole cell lysate

Lane 3 : A549 whole cell lysate

Lane 4 : HepG2 whole cell lysate

Lane 5 : 293T whole cell lysate

Lane 6 : K562 whole cell lysate

Lane 7 : Jurkat whole cell lysate

Loading : 40 μ g

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).