

Calpain S1 antibody [AT1D11]

Cat. No. GTX57626

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, FCM
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Calculated MW 28 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT1D11 is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CAPNS3 protein.
Purification	Protein A Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX57626 FCM Image

FACS analysis of HeLa cells using GTX57626 Calpain S1 antibody.

Cell Number: 1×10^6 cells

Primary antibody: Red line

Antibody amount: 2-5 μ g



GTX57626 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX57626 Calpain S1 antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100


GTX57626 WB Image

WB analysis of various samples using GTX57626 Calpain S1 antibody.

Lane 1 : Recombinant Human CAPNS1 protein (50 ng)

Lane 2 : HeLa whole cell lysate (40 μ g)

Lane 3 : A431 whole cell lysate (40 μ g)

Lane 4 : 293T whole cell lysate (40 μ g)

Lane 5 : Balb/3T3 whole cell lysate (40 μ g)

Lane 6 : U87-MG whole cell lysate (40 μ g)



For full product information, images and publications, please visit our [website](https://www.genetex.com).