

GNMT antibody [AT5D1]

Cat. No. GTX57629

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 33 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT5D1 is derived from hybridization of mouse F2 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human GNMT protein
Purification	Protein A Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX57629 ICC/IF Image

ICC/IF analysis of HepG2 cells using GTX57629 GNMT antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100



GTX57629 WB Image


WB analysis of various samples using GTX57629 GNMT antibody.

Lane 1 : 293T cell lysate

Lane 2 : GNMT transfected 293T cell lysate

Loading : 10 µg

Dilution : 1:1000



GTX57629 WB Image

WB analysis of various samples using GTX57629 GNMT antibody.

Lane 1 : Recombinant Human GNMT protein (50 ng)

Lane 2 : Mouse liver tissue lysate (40 µg)



For full product information, images and publications, please visit our [website](https://www.genetex.com).