

Peroxiredoxin 5 antibody [AT6A10]

Cat. No. GTX57635

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Calculated MW 22 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT6A10 is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PRDX7 protein.
Purification	Protein A Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX57635 ICC/IF Image

ICC/IF analysis of Hep3B cells using GTX57635 Peroxiredoxin 5 antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100




GTX57635 FCM Image

FACS analysis of HeLa cells using GTX57635 Peroxiredoxin 5 antibody.

Cell Number: 1×10^6 cells

Primary antibody: Red line

Antibody amount: 2-5 μ g




GTX57635 ICC/IF Image

ICC/IF analysis of A549 cells using GTX57635 Peroxiredoxin 5 antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100



GTX57635 WB Image

WB analysis of various samples using GTX57635 Peroxiredoxin 5 antibody.

Lane 1 : Recombinant Human PRDX5 protein (50 ng)

Lane 2 : MCF-7 whole cell lysate (40 μ g)

Lane 3 : A549 whole cell lysate (40 μ g)

Lane 4 : HeLa whole cell lysate (40 μ g)

Lane 5 : 293T whole cell lysate (40 μ g)

Lane 6 : Jurkat whole cell lysate (40 μ g)



For full product information, images and publications, please visit our [website](https://www.genetex.com).