

PGAM2 antibody [AT5A7]

Cat. No. GTX57652

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 29 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT5A7 is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PGAM4 protein.
Purification	Protein A Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX57652 ICC/IF Image

ICC/IF analysis of Hep3B cells using GTX57652 PGAM2 antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100



GTX57652 WB Image

WB analysis of various samples using GTX57652 PGAM2 antibody.

Lane 1 : Recombinant Human PGAM2 protein (50 ng)

Lane 2 : HeLa whole cell lysate (40 μ g)

Lane 3 : HepG2 whole cell lysate (40 μ g)

Lane 4 : 293T whole cell lysate (40 μ g)

Lane 5 : Jurkat whole cell lysate (40 μ g)

Lane 6 : NIH-3T3 whole cell lysate (40 μ g)

Lane 7 : A549 whole cell lysate (40 μ g)

Lane 8 : MCF-7 whole cell lysate (40 μ g)



For full product information, images and publications, please visit our [website](https://www.genetex.com).