

GSS antibody [AT12F2]

Cat. No. GTX57671

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT12F2 is derived from hybridization of mouse F2 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human GSS protein
Purification	Protein A Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX57671 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX57671 GSS antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100


GTX57671 WB Image

WB analysis of various samples using GTX57671 GSS antibody.

Lane 1 : MDA-MB231 whole cell lysate

Lane 2 : WiDr whole cell lysate

Lane 3 : SW480 whole cell lysate

Lane 4 : A427 whole cell lysate

Lane 5 : MCF-7 whole cell lysate

Lane 6 : 293T whole cell lysate

Lane 7 : HeLa whole cell lysate

Lane 8 : HT1080 whole cell lysate

Lane 9 : U87-MG whole cell lysate

Lane 10 : A549 whole cell lysate

Loading : 40 µg

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).