

ENO1 + ENO2 antibody [AT1G7]

Cat. No. GTX57673

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, FCM
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Product Note Based on internal testing result, this antibody does not cross react with human ENO3.

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT1G7 is derived from hybridization of mouse F2 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human Alpha-enolase protein.
Purification	Protein A Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX57673 WB Image

WB analysis of various samples using GTX57673 ENO1 + ENO2 antibody [AT1G7].

Lane 1 : PC-3 whole cell lysate


Lane 2 : MCF-7 whole cell lysate

Lane 3 : 293T whole cell lysate

Lane 4 : HeLa whole cell lysate

Loading : 40 µg

Dilution : 1:1000



GTX57673 FCM Image

FACS analysis of LNCap cells using GTX57673 ENO1 + ENO2 antibody [AT1G7].

Cell Number: 1 x 10⁶ cells

Primary antibody: Red line

Antibody amount: 2-5 µg


GTX57673 WB Image

WB analysis of various samples using GTX57673 ENO1 + ENO2 antibody [AT1G7].

Lane 1 : Recombinant Human ENO1 (Alpha-enolase) protein

Lane 2 : Recombinant Human ENO2 (Gamma-enolase) protein

Lane 3 : Recombinant Human ENO3 (Beta-enolase) protein

Loading : 50 ng

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).