

CDC37 antibody [AT3G7]

Cat. No. GTX57677

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 44 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT3G7 is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CDC39 protein.
Purification	Protein A Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX57677 WB Image

WB analysis of various samples using GTX57677 CDC37 antibody.

Lane 1 : Jurkat whole cell lysate
 Lane 2 : A431 whole cell lysate
 Lane 3 : K562 whole cell lysate
 Lane 4 : LNCap whole cell lysate
 Lane 5 : SW480 whole cell lysate
 Lane 6 : HeLa whole cell lysate
 Lane 7 : MCF-7 whole cell lysate
 Lane 8 : 293T whole cell lysate
 Loading : 40 µg
 Dilution : 1:1000


GTX57677 ICC/IF Image

ICC/IF analysis of Balb/3T3 cells using GTX57677 CDC37 antibody.

Blue: DAPI
 Green: Primary antibody
 Dilution: 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).