

nm23-H2 antibody [AT2E3]

Cat. No. GTX57712

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 17 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The clone AT2E3 is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human NME4 protein.
Purification	Protein A Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57712 WB Image

WB analysis of various samples using GTX57712 nm23-H2 antibody.

Lane 1 : HeLa whole cell lysate


Lane 2 : A549 whole cell lysate

Lane 3 : MCF-7 whole cell lysate

Lane 4 : LNCap whole cell lysate

Loading : 40 µg

Dilution : 1:1000



GTX57712 WB Image


WB analysis of various samples using GTX57712 nm23-H2 antibody.

Lane 1 : 293T cell lysate

Lane 2 : NME2 transfected 293T cell lysate

Loading : 5 µg

Dilution : 1:1000



GTX57712 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX57712 nm23-H2 antibody.

Blue: DAPI

Green: Primary antibody

Dilution: 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).