

Hexokinase I + Hexokinase II + Hexokinase III antibody [AT4C12]

Cat. No. GTX57726

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent

Not tested in other applications.

Product Note Based on internal testing result, this antibody does not cross react with human Hexokinase 4.

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human Hexokinase1,2,3 (1-923aa) purified from E. coli.
Purification	Protein A Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX57726 WB Image

WB analysis of various samples using GTX57726 Hexokinase I + Hexokinase II + Hexokinase III antibody [AT4C12].

Lane 1 : HeLa whole cell lysate


Lane 2 : Jurkat whole cell lysate

Lane 3 : K562 whole cell lysate

Lane 4 : Mouse brain tissue lysate

Loading : 40 µg

Dilution : 1:1000



GTX57726 WB Image

WB analysis of various samples using GTX57726 Hexokinase I + Hexokinase II + Hexokinase III antibody [AT4C12].

Lane 1 : Recombinant Human Hexokinase 1 protein

Lane 2 : Recombinant Human Hexokinase 2 protein

Lane 3 : Recombinant Human Hexokinase 3 protein

Lane 4 : Recombinant Human Hexokinase 4 protein

Loading : 20 ng

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).