

EphA1 antibody

Cat. No. GTX57738

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, FCM, LCI
Reactivity	Human, Mouse, Rat

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent
LCI	Assay dependent

Not tested in other applications.

Calculated MW 108 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.8 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide (C)KKEPRQLELTWAGSR, corresponding to amino acid residues 457-471 of human EphA1 (Accession P21709). Extracellular, N-terminus.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTx57738 WB Image

WB analysis of rat (lanes 1 and 3) and mouse (lanes 2 and 4) brain membrane lysates using GTx57738 EphA1 antibody.

Lane 1 and 3 : Primary antibody

Lane 2 and 4 : Primary antibody preincubated with the control peptide antigen

Dilution : 1:200



GTx57738 FCM Image

FACS analysis of Jurkat cells using GTx57738 EphA1 antibody.

Black : Unstained cells

Green : Cells with primary antibody

Dilution : 1:20


GTx57738 LCI Image

Live cell imaging analysis of MCF-7 cells using GTx57738 EphA1 antibody.

Panel A : Primary antibody

Panel B : Live view of the cells

Panel C : Merged of panel A and B

Dilution : 1:25



For full product information, images and publications, please visit our [website](https://www.genetex.com).