

PI3 kinase p110 beta antibody, C-term

Cat. No. GTX59597

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA, IHC
Reactivity	Human

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5-2.5 µg/ml
ELISA	0.25 µg/ml
IHC	5-25 µg/ml

Not tested in other applications.

Calculated MW 122.76 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide selected from the C-terminal region of human PI3Kβ
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx59597 WB Image

WB analysis of Jurkat cell lysate using GTx59597 PI3 kinase p110 beta antibody, C-term.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

Date 2026 / 10 / 09 Page 2 of 2