

Albumin antibody [3F4]

Cat. No. GTX60252

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:300-1000
IHC-P	1:50-400

Not tested in other applications.

Calculated MW 69 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	1% BSA, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human serum albumin
Purification	Protein G purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX60252 IHC-P Image

IHC-P analysis of human hepatocellular carcinoma tissue using GTX60252 Albumin antibody [3F4].

Antigen retrieval : Boiling in sodium citrate buffer (pH6) for 15min

Dilution : 1:500



GTX60252 WB Image

WB analysis of various samples using GTX60252 Albumin antibody [3F4].

Dilution : 1:300

Lane 1 : mouse lung lysates

Lane 2 : mouse spleen lysates


GTX60252 WB Image

WB analysis of A549 cell lysates using GTX60252 Albumin antibody [3F4].

Dilution : 1:300



For full product information, images and publications, please visit our [website](https://www.genetex.com).