

# NGN3 antibody

**Cat. No. GTX60254**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, IHC-P
<b>Reactivity</b>	Human, Mouse, Rat

References ( 1 )

Package

100 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:300-1000
IHC-P	1:50-400

Not tested in other applications.

**Calculated MW** 23 kDa. ( [Note](#) )

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	1% BSA, 50% Glycerol
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human NGN3(60-110).
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Unconjugated

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTx60254 IHC-P Image

IHC-P analysis of rat brain tissue using GTx60254 NGN3 antibody.

Antigen retrieval : 0.4% pepsin at 37°C for 30min


Dilution : 1:200



### GTx60254 WB Image

WB analysis of mouse liver tissue lysate using GTx60254 NGN3 antibody.

Dilution : 1:300




### GTx60254 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTx60254 NGN3 antibody.

Antigen retrieval : Boiling in sodium citrate buffer (pH6) for 15min

Dilution : 1:400



### GTx60254 IHC-P Image

IHC-P analysis of mouse small intestine tissue using GTx60254 NGN3 antibody.

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).