

NGN3 antibody

Cat. No. GTX60254

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

References (1)
 Package
 100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:300-1000
IHC-P	1:50-400

Not tested in other applications.

Calculated MW 23 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	1% BSA, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human NGN3(60-110).
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

**GTX60254 IHC-P Image**

IHC-P analysis of rat brain tissue using GTX60254 NGN3 antibody.


Antigen retrieval : 0.4% pepsin at 37°C for 30min

Dilution : 1:200

**GTX60254 WB Image**

WB analysis of mouse liver tissue lysate using GTX60254 NGN3 antibody.


Dilution : 1:300

**GTX60254 IHC-P Image**

IHC-P analysis of human breast carcinoma tissue using GTX60254 NGN3 antibody.

Antigen retrieval : Boiling in sodium citrate buffer (pH6) for 15min

Dilution : 1:400

**GTX60254 IHC-P Image**

IHC-P analysis of mouse small intestine tissue using GTX60254 NGN3 antibody.

Dilution : 1:200



For full product information, images and publications, please visit our [website](#).