

SEK1 / MKK4 antibody [5H4]

Cat. No. GTX60399

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 44 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human MAP2K4 expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX60399 WB Image

WB analysis of HepG2 (1), K562 (2), and HEK293 (3) cell lysate using GTX60399 SEK1 / MKK4 antibody [5H4].


GTX60399 IHC-P Image


IHC-P analysis of human muscle tissue using GTX60399 SEK1 / MKK4 antibody [5H4].


GTX60399 FCM Image

FACS analysis of K562 cells using GTX60399 SEK1 / MKK4 antibody [5H4].

Green : SEK1 / MKK4

Purple : negative control


GTX60399 ELISA Image

ELISA analysis of antigen using GTX60399 SEK1 / MKK4 antibody [5H4].

Red : Control antigen 100ng

Purple : Antigen 10ng

Green : Antigen 50ng

Blue : Antigen 100ng



For full product information, images and publications, please visit our [website](http://www.genetex.com).