

CD133 antibody [3F10]

Cat. No. GTX60471

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P, FCM, ELISA
Reactivity	Human

References (5)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Synthesized peptide of human CD133.
Purification	Unpurified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX60471 IHC-P Image

IHC-P analysis of human breast cancer tissues (left) and human esophageal cancer tissues (right) using GTX60471 CD133 antibody [3F10]


GTX60471 ELISA Image


ELISA analysis of antigen using GTX60471 CD133 antibody [3F10].

Red : Control antigen 100ng

Purple : Antigen 10ng

Green : Antigen 50ng

Blue : Antigen 100ng


GTX60471 FCM Image

FACS analysis of HeLa cells using GTX60471 CD133 antibody [3F10].

Blue : CD133

Red : negative control



For full product information, images and publications, please visit our [website](http://www.genetex.com).