

Fibrinogen gamma antibody [5A6]

Cat. No. GTX60498

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF, ELISA
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human FGG(AA : 210 - 437) expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX60498 WB Image

WB analysis of human Fibrinogen gamma (AA: 210-437) recombinant protein using GTX60498 Fibrinogen gamma antibody [5A6].



GTX60498 ICC/IF Image

ICC/IF analysis of NIH3T3 cells using GTX60498 Fibrinogen gamma antibody [5A6].

Green : Fibrinogen gamma

Blue: DRAQ5 fluorescent DNA dye

Red: Actin filaments


GTX60498 ELISA Image

ELISA analysis of antigen using GTX60498 Fibrinogen gamma antibody [5A6].

Red : Control antigen 100ng

Purple : Antigen 10ng

Green : Antigen 50ng

Blue : Antigen 100ng



For full product information, images and publications, please visit our [website](http://www.genetex.com).