

## MEKK2 antibody [4B4]

Cat. No. GTX60504

|              |                   |
|--------------|-------------------|
| Host         | Mouse             |
| Clonality    | Monoclonal        |
| Isotype      | IgG1              |
| Applications | WB, ICC/IF, ELISA |
| Reactivity   | Human, Mouse      |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1/500 - 1/2000       |
| ICC/IF             | 1/200 - 1/1000       |
| ELISA              | 1/10000              |

Not tested in other applications.

Calculated MW 70 kDa. ( [Note](#) )

## Properties

|              |                                                                                                                                                                                                                            |
|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form         | Liquid                                                                                                                                                                                                                     |
| Buffer       | Ascites                                                                                                                                                                                                                    |
| Preservative | 0.03% Sodium azide                                                                                                                                                                                                         |
| Storage      | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Immunogen    | Purified recombinant fragment of human MAP3K2 expressed in E. Coli.                                                                                                                                                        |
| Purification | Unpurified                                                                                                                                                                                                                 |
| Conjugation  | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX60504 WB Image**

WB analysis of HEK293 (1) and MEKK2 (AA: 148-359)--hlgGFc transfected HEK293 (2) cell lysate using GTX60504 MEKK2 antibody [4B4].



**GTX60504 ICC/IF Image**

ICC/IF analysis of 3T3-L1 cells using GTX60504 MEKK2 antibody [4B4].

Green : MEKK2

Blue: DRAQ5 fluorescent DNA dye

Red: Actin filaments


**GTX60504 ELISA Image**

ELISA analysis of antigen using GTX60504 MEKK2 antibody [4B4].

Red : Control antigen 100ng

Purple : Antigen 10ng

Green : Antigen 50ng

Blue : Antigen 100ng



For full product information, images and publications, please visit our [website](https://www.genetex.com).