

Alkaline Phosphatase (Tissue Non-Specific) antibody [2F4]

Cat. No. GTX60568

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, ELISA
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 57 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human Alkaline Phosphatase (Tissue Non-Specific) expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX60568 IHC-P Image

IHC-P analysis of ovarian cancer tissue using GTX60568 Alkaline Phosphatase (Tissue Non-Specific) antibody [2F4].




GTX60568 WB Image

WB analysis of HeLa (1), and NTERA-2 (2) cell lysate using GTX60568 Alkaline Phosphatase (Tissue Non-Specific) antibody [2F4].



GTX60568 FCM Image

FACS analysis of MCF-7 cells using GTX60568 Alkaline Phosphatase (Tissue Non-Specific) antibody [2F4].
Green : Alkaline Phosphatase (Tissue Non-Specific)
Red : negative control



GTX60568 ELISA Image

ELISA analysis of antigen using GTX60568 Alkaline Phosphatase (Tissue Non-Specific) antibody [2F4].
Black : Control antigen 100ng
Purple : Antigen 10ng
Blue : Antigen 50ng
Red : Antigen 100ng



For full product information, images and publications, please visit our [website](http://www.genetex.com).