

RAP1A antibody [5F8]

Cat. No. GTX60585

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, FCM, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 21 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human RAP1A expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx60585 FCM Image


FACS analysis of A431 cells using GTx60585 RAP1A antibody [5F8].

Green : RAP1A

Red : negative control


GTx60585 WB Image

WB analysis of HEK293 (1) and RAP1A (AA: 28-180)-hlgGfc transfected HEK293 (2) cell lysate using GTx60585 RAP1A antibody [5F8].


GTx60585 ELISA Image


ELISA analysis of antigen using GTx60585 RAP1A antibody [5F8].

Black : Control antigen 100ng

Purple : Antigen 10ng

Blue : Antigen 50ng

Red : Antigen 100ng


GTx60585 WB Image

WB analysis of human RAP1A (AA: 28-180) recombinant protein using GTx60585 RAP1A antibody [5F8].



For full product information, images and publications, please visit our [website](http://www.genetex.com).