

Hsp90 alpha antibody [5G5]

Cat. No. GTX60587

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	WB, IHC-P, FCM, ELISA
Reactivity	Human, Mouse, Monkey

References (1) Package 100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 85 kDa. (Note)


Properties	
Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human HSP90AA1 expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.


Date 2025 / 11 / 07 Page 1 of 2

DATA IMAGES



GTX60587 IHC-P Image


IHC-P analysis of lung cancer tissue using GTX60587 Hsp90 alpha antibody [5G5].



GTX60587 FCM Image

FACS analysis of HeLa cells using GTX60587 Hsp90 alpha antibody [5G5].

Green: Hsp90 alpha Red: negative control




GTX60587 ELISA Image

ELISA analysis of antigen using GTX60587 Hsp90 alpha antibody [5G5].

Black: Control antigen 100ng

Purple: Antigen 10ng Blue: Antigen 50ng Red: Antigen 100ng



GTX60587 WB Image

WB analysis of human Hsp90 alpha (AA: 275-484) recombinant protein using GTX60587 Hsp90 alpha antibody [5G5].



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 07 Page 2 of 2