

HAS3 antibody [3C9]

Cat. No. GTX60588

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, ELISA
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 63 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human HAS3 expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX60588 WB Image


WB analysis of human HAS3 (AA: 312-364) recombinant protein using GTX60588 HAS3 antibody [3C9].


GTX60588 FCM Image


FACS analysis of HeLa cells using GTX60588 HAS3 antibody [3C9].

Green : HAS3

Red : negative control


GTX60588 IHC-P Image

IHC-P analysis of cervical cancer tissue using GTX60588 HAS3 antibody [3C9].


GTX60588 IHC-P Image

IHC-P analysis of human prostate tissue using GTX60588 HAS3 antibody [3C9].



For full product information, images and publications, please visit our [website](https://www.genetex.com).