

Splunc2 antibody [1F12]

Cat. No. GTX60590

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 25 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of mouse Splunc2 expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX60590 ICC/IF Image

ICC/IF analysis of HepG2 cells using GTX60590 Splunc2 antibody [1F12].


Green : Splunc2

Red: Actin filaments



GTX60590 IHC-P Image

IHC-P analysis of human kidney tissue using GTX60590 Splunc2 antibody [1F12].



GTX60590 FCM Image

FACS analysis of HeLa cells using GTX60590 Splunc2 antibody [1F12].

Green : mouse Splunc2

Red : negative control



For full product information, images and publications, please visit our [website](#).