

PKA C gamma antibody [2E4]

Cat. No. GTX60597

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 40 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.03% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant fragment of human PRKACG expressed in E. Coli.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX60597 WB Image

WB analysis of human PKA C gamma (AA: 164-351) recombinant protein using GTX60597 PKA C gamma antibody [2E4].




GTX60597 FCM Image

FACS analysis of MCF-7 cells using GTX60597 PKA C gamma antibody [2E4].
 Green : PKA C gamma
 Red : negative control



GTX60597 ELISA Image

ELISA analysis of antigen using GTX60597 PKA C gamma antibody [2E4].
 Black : Control antigen 100ng
 Purple : Antigen 10ng
 Blue : Antigen 50ng
 Red : Antigen 100ng



GTX60597 IHC-P Image

IHC-P analysis of rectum cancer tissue using GTX60597 PKA C gamma antibody [2E4].



For full product information, images and publications, please visit our [website](#).