

PI3 kinase p85 alpha antibody [6G10]

Cat. No. GTX60633

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Human, Mouse

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/100
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 84 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human PI3 kinase p85 alpha (AA: 159-388) expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX60633 ICC/IF Image

ICC/IF analysis of HepG2 cells using GTX60633 PI3 kinase p85 alpha antibody [6G10].

Green : PI3 kinase p85 alpha

Blue: DRAQ5 fluorescent DNA dye

Red: Actin filaments




GTX60633 FCM Image

FACS analysis of NIH3T3 cells using GTX60633 PI3 kinase p85 alpha antibody [6G10].

Green : PI3 kinase p85 alpha

Red : negative control



GTX60633 WB Image

WB analysis of human PI3 kinase p85 alpha recombinant protein using GTX60633 PI3 kinase p85 alpha antibody [6G10].



For full product information, images and publications, please visit our [website](#).