

## CSNK2A2 antibody [1E1]

Cat. No. GTX60636

|              |                   |
|--------------|-------------------|
| Host         | Mouse             |
| Clonality    | Monoclonal        |
| Isotype      | IgG1              |
| Applications | WB, FCM, ELISA    |
| Reactivity   | Human, Mouse, Rat |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1/500 - 1/2000       |
| FCM                | 1/200 - 1/400        |
| ELISA              | 1/10000              |

Not tested in other applications.

Calculated MW 41 kDa. ( [Note](#) )

## Properties

|              |                                                                                                                                                                                                                            |
|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form         | Liquid                                                                                                                                                                                                                     |
| Buffer       | Ascites                                                                                                                                                                                                                    |
| Preservative | 0.03% Sodium azide                                                                                                                                                                                                         |
| Storage      | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Immunogen    | Purified recombinant fragment of human CSNK2A2 (AA: 194-350) expressed in E. Coli.                                                                                                                                         |
| Purification | Unpurified                                                                                                                                                                                                                 |
| Conjugation  | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



**GTX60636 WB Image**

WB analysis of HeLa (1), MCF-7 (2), HepG2 (3), Jurkat (4), NIH3T3 (5), and PC-12 (6) cell lysate using GTX60636 CSNK2A2 antibody [1E1].




**GTX60636 FCM Image**

FACS analysis of HepG2 cells using GTX60636 CSNK2A2 antibody [1E1].

Green : CSNK2A2

Red : negative control



**GTX60636 ELISA Image**

ELISA analysis of antigen using GTX60636 CSNK2A2 antibody [1E1].

Black : Control antigen 100ng

Purple : Antigen 10ng

Blue : Antigen 50ng

Red : Antigen 100ng



For full product information, images and publications, please visit our [website](https://www.genetex.com).