

ID2 antibody [4E12G5]

Cat. No. GTX60689

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 15 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human ID2 (AA: 1-134) expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX60689 IHC-P Image

IHC-P analysis of rectum cancer tissue using GTX60689 ID2 antibody [4E12G5].



GTX60689 IHC-P Image

IHC-P analysis of breast cancer tissue using GTX60689 ID2 antibody [4E12G5].




GTX60689 FCM Image

FACS analysis of SK-N-SH cells using GTX60689 ID2 antibody [4E12G5].

Green : ID2

Purple : negative control



GTX60689 WB Image

WB analysis of human ID2 recombinant protein using GTX60689 ID2 antibody [4E12G5].



For full product information, images and publications, please visit our [website](https://www.genetex.com).