

SLC27A5 antibody [4B11C10]

Cat. No. GTX60693

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Calculated MW 75 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human SLC27A5 (AA: 508-570) expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX60693 IHC-P Image

IHC-P analysis of human stomach tissue using GTX60693 SLC27A5 antibody [4B11C10].



GTX60693 ELISA Image


ELISA analysis of antigen using GTX60693 SLC27A5 antibody [4B11C10].

Black : Control antigen 100ng

Purple : Antigen 10ng


Blue : Antigen 50ng

Red : Antigen 100ng



GTX60693 IHC-P Image

IHC-P analysis of liver cancer tissue using GTX60693 SLC27A5 antibody [4B11C10].



GTX60693 WB Image

WB analysis of 3T3L1 (1) and NIH3T3 (2) cell lysate using GTX60693 SLC27A5 antibody [4B11C10].



For full product information, images and publications, please visit our [website](http://www.genetex.com).