

WHSC2 antibody [6B11H8]

Cat. No. GTX60733

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Human, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 57 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human WHSC2 (AA: 280-511) expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx60733 IHC-P Image

IHC-P analysis of colon cancer tissue using GTX60733 WHSC2 antibody [6B11H8].



GTx60733 ELISA Image


ELISA analysis of antigen using GTX60733 WHSC2 antibody [6B11H8].

Black : Control antigen 100ng

Purple : Antigen 10ng


Blue : Antigen 50ng

Red : Antigen 100ng



GTx60733 WB Image

WB analysis of Jurkat (1), HeLa (2), HEK293 (3), A549 (5), SPC-A-1 (6) cell lysate, and Rat brain (4) tissue lysate using GTX60733 WHSC2 antibody [6B11H8].



GTx60733 IHC-P Image

IHC-P analysis of ovarian cancer tissue using GTX60733 WHSC2 antibody [6B11H8].



For full product information, images and publications, please visit our [website](http://www.genetex.com).