

PHC1 antibody [1F3F3]

Cat. No. GTX60752

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, FCM, ELISA
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 106 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human PHC1 (AA: 758-1004) expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX60752 ELISA Image


ELISA analysis of antigen using GTX60752 PHC1 antibody [1F3F3].

Black : Control antigen 100ng

Purple : Antigen 10ng

Blue : Antigen 50ng


Red : Antigen 100ng


GTX60752 FCM Image

FACS analysis of HEK293 cells using GTX60752 PHC1 antibody [1F3F3].

Green : PHC1

Red : negative control


GTX60752 WB Image

WB analysis of HEK293 (1) and PHC1 (AA: 758-1004)-hIgGFc transfected HEK293 (2) cell lysate using GTX60752 PHC1 antibody [1F3F3].



For full product information, images and publications, please visit our [website](http://www.genetex.com).