

GATA6 antibody [2F10G3]

Cat. No. GTX60772

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, FCM, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/50 - 1/100
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 60 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human GATA6 (AA: 491-557) expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX60772 ELISA Image


ELISA analysis of antigen using GTX60772 GATA6 antibody [2F10G3].

Black : Control antigen 100ng


Purple : Antigen 10ng

Blue : Antigen 50ng

Red : Antigen 100ng


GTX60772 WB Image

WB analysis of HEK293 (1) and GATA6 (AA: 491-557)-hlgGfc transfected HEK293 (2) cell lysate using GTX60772 GATA6 antibody [2F10G3].


GTX60772 FCM Image

FACS analysis of HEK293 cells using GTX60772 GATA6 antibody [2F10G3].

Green : GATA6

Red : negative control



For full product information, images and publications, please visit our [website](http://www.genetex.com).