

BRCA1 antibody [6C6D2]

Cat. No. GTX60779

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
IHC-P	1/200 - 1/1000
ELISA	1/10000

Not tested in other applications.

Calculated MW 208 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human BRCA1 (AA: 229-335) expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX60779 IHC-P Image

IHC-P analysis of human esophagus tissue using GTX60779 BRCA1 antibody [6C6D2].




GTX60779 IHC-P Image

IHC-P analysis of esophagus cancer tissue using GTX60779 BRCA1 antibody [6C6D2].



GTX60779 WB Image

WB analysis of HEK293 (1) and BRCA1 (AA: 229-335)-hlgGfC transfected HEK293 (2) cell lysate using GTX60779 BRCA1 antibody [6C6D2].



GTX60779 ELISA Image

ELISA analysis of antigen using GTX60779 BRCA1 antibody [6C6D2].

Black : Control antigen 100ng

Purple : Antigen 10ng

Blue : Antigen 50ng

Red : Antigen 100ng



For full product information, images and publications, please visit our [website](http://www.genetex.com).