

IL3 Receptor alpha antibody [8E11C5]

Cat. No. GTX60786

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, FCM, ELISA
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1/500 - 1/2000
ICC/IF	1/200 - 1/1000
IHC-P	1/200 - 1/1000
FCM	1/200 - 1/400
ELISA	1/10000

Not tested in other applications.

Calculated MW 43 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human IL3RA (AA: 200-305) expressed in E. Coli.
Purification	Protein G Purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX60786 WB Image

WB analysis of HEK293 (1) and IL3RA (AA: 200-305)-hlgGFc transfected HEK293 (2) cell lysate using GTX60786 IL3 Receptor alpha antibody [8E11C5].



GTX60786 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX60786 IL3 Receptor alpha antibody [8E11C5].

Green : IL3 Receptor alpha

Blue: DRAQ5 fluorescent DNA dye

Red: Actin filaments


GTX60786 FCM Image

FACS analysis of HeLa cells using GTX60786 IL3 Receptor alpha antibody [8E11C5].

Green : IL3 Receptor alpha

Red : negative control



For full product information, images and publications, please visit our [website](https://www.genetex.com).