

VWF antibody

Cat. No. GTX60934

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Mouse, Rat

References (3)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.5µg/ml
IHC-P	0.5-1µg/ml

Not tested in other applications.

Calculated MW 309 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1% Na ₂ HPO ₄ , 0.45% NaCl, 2.5% BSA
Preservative	0.025% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	500 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	E.coli-derived mouse VWF recombinant protein (Position: M1304-E1452).
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX60934 IHC-P Image

IHC-P analysis of mouse liver tissue using GTX60934 VWF antibody.


GTX60934 WB Image

WB analysis of 50ug mouse lung tissue lysate using GTX60934 VWF antibody.


Dilution : 0.5 µg/mL


GTX60934 IHC-P Image

IHC-P analysis of rat lung tissue using GTX60934 VWF antibody.

Antigen retrieval : Heat mediated antigen retrieval was performed in citrate buffer (pH6, epitope retrieval solution) for 20 mins

Dilution : 1µg/ml


GTX60934 WB Image

WB analysis of 0.5ng recombinant mouse VWF protein using GTX60934 VWF antibody.

Dilution : 0.5 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).