

GLP1 antibody [03]

Cat. No. GTX60954

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 21 kDa. ([Note](#))

Product Note This antibody binds GLP-1(1-37), GLP-1(7-37) and GLP-1(9-37).Cross-reaction with C-terminally amidated GLP-1 is 0-16% in various circumstances.

Properties

Form	Liquid
Buffer	10mM Phosphate, 500mM NaCl
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic GLP-1(9-37)
Purification	Protein-A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX60954 ELISA Image

The calibration curve of a sandwich assay for Non-amidated GLP-1 using GTX60954 as the capture antibody and GTX60952 -02 as the biotinylated detection antibody.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 2 of 2