

## Hemopexin antibody [4] (Biotin)

Cat. No. GTX60959-02

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ELISA
Reactivity	Human

Package  
50 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

**Product Note** This antibody is specific for human hemopexin.

## Properties

Form	Liquid
Buffer	10mM Phosphate, 140mM NaCl
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Hemopexin purified from human plasma
Purification	Biotinylated
Conjugation	Biotin

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES

**GTX60959-02 ELISA Image**

The calibration curve of a sandwich assay for Hemopexin using GTX60981 as capture antibody and GTX60959 -02 as the biotinylated detection antibody.



For full product information, images and publications, please visit our [website](#).