

PAPP A antibody [9]

Cat. No. GTX60974

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IHC-P, IHC-Fr, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	181 kDa. (Note)
Product Note	This antibody is specific for native human PAPP-A. Heat and SDS denatured PAPP-A is also recognized, but not acid treated PAPP-A.

Properties

Form	Liquid
Buffer	10mM Phosphate, 500mM NaCl
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified human PAPP-A
Purification	Protein-A purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX60974 IHC-Fr Image

Staining is restricted to the syncytiotrophoblast of the chorionic villi in frozen human placenta tissue, stained with GTX60974.



GTX60974 IHC-P Image

Staining is restricted to the syncytiotrophoblast of the chorionic villi in paraffin embedded human placenta tissue, stained with GTX60974.



For full product information, images and publications, please visit our [website](#).