

Leptin antibody [11] (Biotin)

Cat. No. GTX60983-02

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA
Reactivity	Human

Package
50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Product Note This antibody binds human leptin.

Properties

Form	Liquid
Buffer	10mM Phosphate, 140mM NaCl
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human leptin
Purification	Biotinylated
Conjugation	Biotin

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX60983-02 ELISA Image**

The calibration curve of a sandwich assay for Leptin using GTX60964 as the capture antibody and GTX60983-02 as the biotinylated detection antibody.



For full product information, images and publications, please visit our [website](#).