

MEK6 antibody [EP558Y], C-term

Cat. No. GTX61432

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at a dilution of 1:10,000

For IHC: Use at a dilution of 1:100 - 250.

Optimal working dilution for a specific application should be determined by the investigator.


Calculated MW 37 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	50mM Tris-Glycine, 150mM NaCl, 0.05% BSA, 40% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	A synthetic peptide corresponding to residues in the C-term of human MEK6 was used as immunogen.
Purification	Tissue culture supernatant
Conjugation	Unconjugated
Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. RabMAb® technology is covered by the following U.S. Patents, No. 5,675,063 and/or 7,429,488.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX61432 WB Image

A. Western blot analysis on Hela cell lysate using anti-MKK6 RabMAb (cat. # GTX61432) dilution 1:10000.



GTX61432 IHC-P Image

B. Immunohistochemical analysis of paraffin-embedded human heart muscle using anti-MKK6 RabMAb (cat. # GTX61432).



For full product information, images and publications, please visit our [website](https://www.genetex.com).